

# RT-PCR enzyme/master mix

- Actively testing enzyme capabilities for one-step RT-qPCR detection of SARS-CoV-2:
  - RapiDxFire qPCR 5X Master Mix GF, ([Catalog # 30050-1, 1mL](#))
  - EpiScript RNase H- Reverse Transcriptase ([Catalog # ERT12925K, 25,000 U](#))
    - Sensitive PCR enzymes & mixes manufactured in ISO13485 qualified facility
  - 2019-nCoV CDC Probe and Primer Kit for SARS-CoV-2 (Catalog # [KIT-nCoV-PP1-1000](#))
  - AccuPlex<sup>™</sup> Coronavirus SARS-CoV-2 Reference Material ([Catalog # 0505-0126](#))
  - Preliminary results look promising

